

turni cek

by Martha Mozartha

Submission date: 19-Jun-2024 10:48PM (UTC+0700)

Submission ID: 2399515669

File name: ABMC_INSISIVA_MM_2024.docx (112.87K)

Word count: 2151

Character count: 13272



Research Article

Hardness comparison of single-shade composite vs multi-shade dental resin nano-composite

Ananda¹, Laura Nuradnin¹, Martha Mozartha², Bebby Arisya Kesumaputri², Danica Anastasia³

¹Dentistry Study Program, Faculty of Medicine, Universitas Sriwijaya, Palembang, 30125, Indonesia

²Department of Biomaterial, Dentistry Study Program, Faculty of Medicine, Universitas Sriwijaya, Palembang, 30125, Indonesia

³Department of Conservative Dentistry, Dentistry Study Program, Faculty of Medicine, Universitas Sriwijaya, Palembang, 30125, Indonesia

Received date: ...; revised date: ...; accepted: ...

DOI:

Abstract

Background: Shade matching of multi-shade dental resin nano-composite is a process that requires time and operator expertise. Single-shade universal composite resin claimed to match all tooth colors to simplify the process. **Objective:** to evaluate the surface hardness of single-shade universal composite compared multi-shade dental resin nano-composite. **Method:** This experimental study tested two composite restorative materials: single-shade universal composite (Omnichroma, Tokuyama Dental, Japan), and nanofilled composite resin (Filtek™ Z350 XT, 3M/ESPE, USA). Thirty specimens divided into 2 groups (n=15) were fabricated by placing the composites into a stainless steel mold (2 mm thick × 12 mm diameter), and photoactivated for 20 seconds, followed by finishing and polishing. Surface hardness was evaluated using a Vickers hardness testing machine using diamond indenter with a load of 1000 g applied for 15s. Data were analyzed by using an independent t-test. **Results:** The mean (standard deviation) surface hardness of single-shade universal and multi-shade dental resin nano-composite were 21.04±4.92 and 19.33±5.61 VHN, respectively. Single-shade universal composite showed no significant difference in surface hardness to multi-shade dental resin nano-composite (p-value > 0.05). **Conclusion:** The surface hardness value of single-shade universal composite resin was comparable to multi-shade dental resin nano-composite.

Keywords: multi-shade dental resin nano-composite; surface hardness; single-shade universal composite; Vickers hardness

INTRODUCTION

The development of restorative materials has progressed significantly, and composite resins are increasingly favored because of their continuous improvements and expanding range of applications. These resins are categorized as hybrid, microhybrid, and microfilled, depending on their filler particle size. The emergence of nanotechnology has paved the way for nanofilled composites, which are now utilized in clinical settings. There are two types of composites with nano-scale fillers: nanofilled and nanohybrid. Nanofilled composites consist entirely of nanometer-sized particles, while nanohybrid composites are a mixture of nanofiller particles and larger-sized particles[1].

The main goal of adding nanosized fillers to resin composite is to enhance the mechanical characteristics of the materials, thus making them suitable for use in both anterior and posterior teeth, and have adequate mechanical properties to function in high-stress-bearing areas. Various laboratory tests have previously been conducted to assess the longevity of restorations in the oral environment, including surface hardness testing [2]. Previous studies showed that the nanocomposites had significantly higher hardness values compared to the hybrid and microfilled composites that were tested due to higher filler content, denser filler particles, and resin content [3,4].

In addition to their excellent mechanical properties, dental nanocomposites also offer advantages such as good aesthetics, low surface roughness, and ease of polishing, resulting in a high-gloss finish due to their small particle size.[5]. Nevertheless, achieving accurate shade matching remains a significant challenge even for highly skilled operators. To avoid this issue, a layering technique is used with several resin composite shades. This results in a wide selection of shades that the clinician must keep in stock.

Saegusa et al. assessed the color matching of structural colored resin composites in comparison to conventional multi shade nano-filled resin composites that contain pigments. The new restorative material, which incorporates smart chromatic technology, has recently been introduced to the market and is claimed to be capable of matching all tooth colors. The findings reveal that the structural colored resin composite outperforms traditional resin composites in color matching ability.[6] This is in line with Cruz da Silva et al., who found that the structural colored resin composite, also known as the universal single-shade composite (Omnichroma), demonstrates better color matching compared to multishade resins.[7] Previous studies

suggest that the use of this single-shade universal composite simplifies the shade matching procedure and shortens treatment time.[8]

Comprehensive testing is necessary to determine the characteristics and clinical performance of this relative new single-shade universal composite. Further studies are needed to evaluate surface hardness in comparison with nanofilled resin composite. The current study was conducted to evaluate the surface hardness of a single-shade composite compared to multi-shade dental resin nano-composite.

MATERIAL AND METHOD

This in vitro study with the post-test-only control group method was conducted at the Dental Hospital of South Sumatera Province and the Mechanical Engineering Laboratory Universitas Sriwijaya in Palembang, Indonesia. Two different resin based composites were investigated: a nanofilled (Filtek™ Z-350 XT, 3M Dental Products, St Paul, MN, USA) and a universal single shade composite (Omnichroma®, Tokuyama Dental, Tokyo, Japan). Table 1 provides an overview of the restorative materials examined in the current study.

Table 1. Products information for tested dental composites[9,10]

Brand	Manufacturer	Type	Resin matrix	Fillers
Filtek™ Z-350 XT	3M ESPE (St. Paul, MN, USA).	Multi-shade dental resin nano-composite	Bis-GMA, Bis-EMA, UDMA, TEGDMA, PEGDMA	20 nm nanosilica fillers, 5.00–20.00 nm agglomerates zirconia/silica particles, 0.60–1.40 µm clusters particle size (78.5% by weight)
Omnichroma®	Tokuyama Dental, Tokyo, Japan)	Single-shade universal composite	UDMA, TEGDMA, mequinol, dibutylhidroxitolueno e UV absorber	79% by weight (68% by volume) of spherical zirconium silica filler (average particle size 0.3 µm, range 0.2–0.4 µm)

Brand	Manufacturer	Type	Organic matrix	Fillers
Filtek™ Z-350 XT	3M ESPE (St. Paul, MN, USA).	Nanofilled universal resin composite	Bis-GMA, Bis-EMA, UDMA, TEGDMA, PEGDMA	20 nm nanosilica fillers, 5.00-20.00 nm agglomerates zirconia/silica particles, 0.60-1.40 µm clusters particle size (78.5% by weight)
Omnichroma®	Tokuyama Dental, Tokyo, Japan)	Single-shade universal composite	UDMA, TEGDMA, mequinol, dibutylhidroxitoluene e UV absorber	79% by weight (68% by volume) of spherical zirconium silica filler (average particle size 0.3 µm, range 0.2-0.4 µm)

The stainless steel mold was placed on a mylar strip set on a glass slide. Each sample was created by placing the composite inside the mold and covering it with another mylar strip. A glass sheet was then used to apply pressure, compressing the composite resin and removing any excess material. The LED Light Curing Unit (LUX I, DTE®) was then used to irradiate each specimen from both the top and bottom for 20 seconds without removing the glass sheet, then the cured specimens were dislodged.

After finishing and polishing using polishing kits, each sample was placed in a small container and labeled with its respective group: Group 1 for nanofilled composite and Group 2 for universal single-shade composite (with n=15 samples in each group). Vickers hardness values were measured using a Vickers hardness tester

machine (Albert Gnehm, Switzerland) with a 1000 g load applied using a diamond tip at the center of the material for a duration of 15 seconds. A microscope was used to measure the diagonals of the square traces. The collected data were statistically analyzed to assess their hardness values.

RESULT

The result of Kolmogorov-Smirnov test indicating a normal distribution ($p > 0.05$). Assessment of any significant differences between the tested composites was performed using the Independent t-test. The measurement results for the surface hardness of these two types of resin composites are presented in Table 2.

Table 2. Surface Hardness Value

Property	Mean±SD		p-value
	Nanofilled	Single-Shade Universal	
Hardness	19,33 ± 5,61	21,04 ± 4,92	0,384

Table 2 shows that universal single shade resin composite has greater surface hardness compared to conventional nanofilled resin composite, but the difference was not statistically significant ($p > 0.05$).

DISCUSSION

Surface hardness testing is advantageous to evaluate a material's properties, including wear resistance to indentation of another harder material. It stands as one of the crucial factors that determine the long-term clinical

performance of resin composite restorative materials when functioned in the oral environment. By measuring surface hardness, one can determine the curing depth of composite resin materials. Rezaei et al. discussed that the surface hardness of composite resin can be correlated with many factors, including the size, type, and distribution of the fillers, the type and viscosity of monomers, the initiator system, and the light-curing unit device [11].

In this study, two types of composite resins were tested: FiltekTM Z 350 XT, a multi-shade nanofilled composite, and Omnicroma, a universal single-shade composite. Based on the results, the two composites exhibited similar surface hardness values. Independent t-tests showed no statistically significant difference between the two. It can be attributed to the similarity in filler particle characteristics between the two composites. Omnicroma and FiltekTM Z 350 XT has a filler loading of 79 wt% and 78.5 wt%, respectively[12].

Kundie et al. reviewed articles that investigated the influence of the fillers size on mechanical properties of composite resins. The findings revealed that filler size plays a role in determining the hardness of composites. Dental composites containing nanoparticles exhibited elevated hardness levels.[13] Omnicroma, which was examined, contains uniformly sized supra-nano spherical filler particles measuring 260 nm in diameter. Similarly, Filtek Z 350 XT contains nanometer-scale spherical particles in the form of clusters and agglomerates. Furthermore, both of them also contain silica and zirconia filler particles [9,12]. Silica fillers are the most widely used in commercial resin fillings, and zirconium oxide is added to further improve the mechanical properties, including hardness. Previous studies have emphasized that incorporating zirconia as a filler in composite materials enhances both strength and toughness, but a decrease occurs when zirconia-silica concentrations are increased. Currently, it is unknown to what extent zirconia-silica concentrations are present in Omnicroma.[14]

The measurement of the initial Vickers hardness values of the bottom and the top surfaces of the composites can provide insight into the depth of polymerization. A Bottom/Top average ratio approaching 1 indicates adequate light transmission from top to bottom and ensured that the base has adequately been polymerized. However, there is a limitation in this study, where

Vickers testing were only conducted on the top surface. Previous research suggests that Vickers hardness values of composites decreased in all the samples after undergoing thermocycling to simulate the oral environment.[15] Further research¹³ is required to determine the aging effects on the surface hardness of Omnicroma.

CONCLUSION

The universal single-shade composite resin (Omnicroma®) exhibits surface hardness values comparable to that of the multi-shade dental resin nano-composite (FiltekTM Z350 XT).

16

CONFLICT OF INTEREST

The authors declare no potential conflict of interest.

REFERENCE

- [1]. Alzraikat H, Burrow M, Maghaireh G, Taha N. Nanofilled resin composite properties and clinical performance: A Review. Operative Dentistry. 2018;43(4). doi:10.2341/17-208-t
- [2]. Nithya K, Sridevi K, Keerthi V, Ravishankar, P. Evaluation of surface roughness, hardness, and gloss of composites after three different finishing and polishing techniques: An in vitro study. Cureus. 2020; doi:10.7759/cureus.7037
- [3]. Abuelenain DA, Neel EAA, Al-Dharab A. Surface and Mechanical Properties of Different Dental Composites. Austin J Dent. 2015;2(2): 1019
- [4]. Poggio C, Lombardini M, Gaviati S, Chiesa M. Evaluation of Vickers hardness and depth of cure of six composite resins photo-activated with different polymerization modes. Journal of Conservative Dentistry. 2012;15(3)237-241
- [5]. Rai R, Gupta R. In vitro evaluation of the effect of two finishing and polishing systems on four esthetic restorative materials. J Conserv Dent. 2013;16:564-7
- [6]. Saegusa M, Kurokawa H, Takahashi N, Takamizawa T, Ishii R, Shiratsuchi K, Miyazaki M. Evaluation of Color-matching Ability of a Structural Colored Resin Composite. Oper Dent. 2021 May 1;46(3):306-315. doi: 10.2341/20-002-L. PMID: 34411249.
- [7]. Cruz da Silva ET, Charamba Leal CF, Miranda SB, Evangelista Santos M, Saeger Meireles S, Maciel de Andrade AK, Japiassú Resende Montes MA. Evaluation of Single-Shade Composite Resin Color Matching on Extracted Human Teeth. Scientific World Journal. 2023 Jun 26;4376545. doi: 10.1155/2023/4376545.

- [8]. Sharma N, Shamant PS. Omnichroma: The See-It -To-Believe -It Technology. *EAS J Dent Oral Med*. 2021 May-Jun; Vol-3, Iss- 3: 100-104. doi:10.36349/easjdom.2021.v03i03.008
- [9]. Furusawa K, Kobayashi S, Yamashita A, Tichy A, Hosaka K, Shimada Y, et al. Effect of filler load on structural coloration and color adjustment potential of resin composites. *Dental Materials Journal*. 2023;42(3):343–50. doi:10.4012/dmj.2022-199
- [10]. Jassá FF, Braga Borges CH, Tonetto MR, de Souza Rastelli AN, Bagnato VS, de Campos EA, Soares Dos Santos RS, Bandeca MC, Cury Saad JR. Long-term surface hardness and monomer conversion of a nanofilled and a microhybrid composite resin. *J Contemp Dent Pract*. 2013 Sep 1;14(5):876-82. doi: 10.5005/jp-journals-10024-1419. PMID: 24685791.
- [11]. Rezaei S, Abbasi M, Sadeghi Mahounak F, Moradi Z. Curing depth and degree of conversion of five bulk-fill composite resins compared to a conventional composite. *The Open Dentistry Journal*. 2019;13(1):422–9. doi:10.2174/1874210601913010422
- [12]. 9. Abreu, et al. Analysis of the color matching of universal resin composites in anterior restorations. *J Esthet Restor Dent*. 2020: 1-8.
- [13]. Kundie, F. et al. (2018). Effects of filler size on the mechanical properties of polymer-filled dental composites: A review of recent developments. *J. Phys. Sci.*, 29(1), 141–165, <https://doi.org/10.21315/jps2018.29.1.10To>
- [14]. Dai S, Chen Y, Yang J, He F, Chen C, Xie H. Surface Treatment Of Nanozirconia Fillers To Strengthen Dental Bisphenol A-Glycidyl Methacrylate-Based Resin Composites. *Int J Nanomedicine*. 2019 Nov 26;14:9185-9197. doi: 10.2147/IJN.S22339
- [15]. Gönder HY, Fidan M. Effect of different polymerization times on color change, translucency parameter, and surface hardness of bulk-fill resin composites. *Niger J Clin Pract*. 2022;25:1751-7

turni cek

ORIGINALITY REPORT

17%

SIMILARITY INDEX

15%

INTERNET SOURCES

14%

PUBLICATIONS

4%

STUDENT PAPERS

PRIMARY SOURCES

1	journalcra.com Internet Source	2%
2	www.hindawi.com Internet Source	2%
3	Submitted to Universitas Muhammadiyah Surakarta Student Paper	2%
4	iopscience.iop.org Internet Source	2%
5	bmcoralhealth.biomedcentral.com Internet Source	2%
6	dentjods.sums.ac.ir Internet Source	1%
7	"Full Issue PDF", Operative Dentistry, 2019 Publication	1%
8	"Abstracts of the 10th Virtual Conseuro 2021 Congress", Clinical Oral Investigations, 2021 Publication	1%
9	smujo.id	

10

dergipark.org.tr

Internet Source

1 %

11

jcd.org.in

Internet Source

1 %

12

Fatma Aytaç Bal, Merve Ağaccioğlu, Osman Demir. "A comparative study of the 5 mm-layer Vickers hardness model with bulk-fill resin-based composites", International Dental Research, 2021

Publication

1 %

13

Loulwa M. Al-Saud, Lina M. Alolyet, Dhayah S. Alenezi. "The Effects of Selected Mouthwashes on the Surface Microhardness of a Single-shade Universal Resin Composite: In Vitro Study", Journal of Advanced Oral Research, 2022

Publication

1 %

14

Maria Teresa Xavier, Ana Luísa Costa, Francisco José Caramelo, Paulo Jorge Palma, João Carlos Ramos. "Evaluation of the Interfaces between Restorative and Regenerative Biomaterials Used in Vital Pulp Therapy", Materials, 2021

Publication

1 %

15

S. H. Dickens, J. W. Stansbury, K. M. Choi, C. J. E. Floyd. "Photopolymerization Kinetics of Methacrylate Dental Resins",
Macromolecules, 2003

Publication

1 %

16

pubmed.ncbi.nlm.nih.gov

Internet Source

1 %

Exclude quotes On

Exclude matches < 1%

Exclude bibliography On